

# EluVac™ Vacuum Manifold for Solid Phase Extraction

The EluVac™ Vacuum Manifold allows simultaneously processing of up to 20 samples under vacuum. It is the ideal solution for enriching, concentrating and isolating analytes in the solid phase extraction.

It is mechanically and chemically robust, flexible and can optimally be used for a wide range of matrices in the environmental, food and feed analysis as well as in pharmaceutical and clinical applications.

A special highlight: By simply turning the lid you can switch from “waste” to “collect” position without disassembling the device.



### Further advantages at a glance:

- + Rack for large volume reservoirs
- + High-quality and resistant materials, allowing the use of most solvents
- + Universally applicable, can easily be adapted to a multitude of purposes
- + Also suitable for 100mL rotary evaporator bottles
- + Saves time and money



For further information  
contact [mail@arcsciences.com](mailto:mail@arcsciences.com)

**The right equipment for special applications:**



Sample rack for 10x 60mL bottles



Sample rack for 20x 4mL vials



Sample rack for 6x 100mL round bottom bottles for rotary evaporator (NS 29/32)



**The ideal completion:**

**Chemically resistant vacuum diaphragm pump, one-stage, type ME 2C**

The modern alternative to water jet pumps for evacuation

- + Continuous, oil-free pumping of gases
- + Chemically resistant materials
- + Very low noise level
- + Compact design
- + Long lifetime
- + Easy-to-change diaphragms and valves
- + Ready-to-operate, incl. 2m PVC-hose

## Ordering information

**EluVac™, Vacuum Manifold, 20 positions**  
Part No. 11148

- Rack for up to 10 sample reservoirs  
Part No. 11098
- Collection rack for 60 mL bottles,  
10 positions  
Part No. 11118
- Collection rack for pear shaped rotary  
evaporator bottles, 100 mL, 6 positions  
Part No. 11112
- Collection rack for 4mL vials, 20 positions  
Part No. 11105

*Further racks on request*

**Complete sets for usage with LCTech  
Elufix™ columns:**

**EluVac™, Vacuum Manifold Set, 20 positions**  
EluVac™ Vacuum Manifold, sample rack and  
sample rack for 60 mL bottles  
Part No. 11146

**EluVac™, Vacuum Manifold Set, 20 positions**  
EluVac™ Vacuum Manifold, sample rack and  
sample rack for rotary evaporator bottles  
Part No. 11147

**ME 2C, vacuum diaphragm pump,  
one-stage**  
230 V, 50-60 Hz, CEE cable  
Part No. 11125

*Other voltages on request*